



Hillrom™



## Hillrom Surgical Tables

Putting time back in your hands



Uniform coupling system across the range



Long service intervals – less downtime, more revenue



Unprecedented imaging capability in mobile tables



Combining technologies for efficiency and cost optimisation



One touch position recall

Hillrom is a part of Baxter

**Baxter**




## WITH A STRONG FOUNDATION, ANYTHING IS POSSIBLE

Our versatile, easy-to-use table platform sustains the volume and variety of procedures you face every day, by helping you support your patients and streamline your workflow. Because the right surgical table isn't just a table – it's the foundation on which care happens.



### Intuitive – quick set up

Hillrom Surgical Tables are quick to set up and easy to use with lightweight components, self-aligning hook coupler system, and an intuitive remote control.



### 360° Imaging simplified

Hillrom Surgical Tables with radiolucent tabletop and carbon fiber equipment make 360° imaging simple. You can also position a C-arm where needed with convenient longitudinal sliding range.



### Multi-functional – eliminating need for dedicated tables

No more need for specialised, single-function tables. Hillrom Surgical Tables easily accommodate vascular, spine and neuro procedures with their flexible configurations and accessories.



### Innovative Light Messaging System

Colour coded LED lights communicate table status – such as brake settings, battery levels and error states, helping to reduce alarm fatigue in your clinical staff.

The best surgical care begins with a strong foundation of safety and support. To learn more call Hillrom Customer Care on 1800 445 576 or visit [hillrom.com.au](http://hillrom.com.au)

## ALWAYS READ THE LABEL AND FOLLOW THE DIRECTIONS FOR USE

CONTACT YOUR HILLROM REPRESENTATIVE FOR THE LATEST VERSIONS

AUTHORISED AUSTRALIAN SPONSOR Hill-Rom Pty. Ltd. Suite 4.01, 2 Lyonpark Road, Macquarie Park NSW Australia 2113

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale of its products.