



**Baxter**

**Welch Allyn**  
**Spot Vital Signs**  
MONITOR

**SIMPLICITY,  
DELIVERED**

The **Spot Vital Signs 4400** device offers an efficient way to capture, access and document patient vitals so you can focus on what really matters — spending more time with your patients.

---

## SIMPLE

An intuitive, all-in-one design makes the **Spot Vital Signs 4400** device easy-to-use for caregivers at any level. Add it to your existing workflow without disruption and begin exams quickly from the touchscreen.

## ACCURATE

Trusted for clinical accuracy and consistent measurements,<sup>1</sup> the **Spot Vital Signs 4400** device can help you make informed treatment decisions at the point of care. The portable design can also help

standardise vitals capture throughout your facility and increase productivity for high-volume environments. Plus, blood pressure averaging is now part of the main workflow, enabling you to take and average multiple readings for improved hypertension detection.<sup>2</sup>

## CONNECTED

Save time and help reduce data entry errors<sup>3</sup> by sending patient data directly to the EMR.



EASY TO NAVIGATE,  
EASY-TO-USE



BP AVERAGING  
WITH ONE TAP



ALL-IN-ONE  
VITALS CAPTURE



EMR CONNECTIVITY  
VIA USB CABLE

# YOUR NEW SOLUTION FOR DETECTING HYPERTENSION

---

## EXPERIENCE BLOOD PRESSURE READINGS YOU CAN BELIEVE IN.

You may have heard about blood pressure averaging, but do you know how it can help improve hypertension screening? Relying on a single blood pressure measurement can be misleading.<sup>2</sup>

First-time readings can yield higher results than subsequent blood pressure measurements.<sup>2</sup> Additionally, it's common

for variations to occur between readings captured in a series.<sup>4</sup> That's why the American Heart Association recommends automated office blood pressure solutions that can capture and average multiple readings during a single visit.<sup>2</sup>

Blood pressure averaging helps provide a more accurate view of your patient's heart health, so you can detect hypertension with confidence.<sup>2</sup>



## ORDERING INFORMATION\*

Part Number	Description
44XT-6	<b>Spot</b> monitor with NIBP and <b>SureTemp</b> Thermometer
44WT-6	<b>Spot</b> monitor with NIBP, <b>SureTemp</b> Thermometer and <b>Nonin</b> SpO <sub>2</sub>
4400-DST	<b>Spot</b> monitor Desk Stand
4400-MBS	<b>Spot</b> monitor Mobile Stand
777-PM2WCX-6	Wall System 77716, <b>PanOptic</b> Ophthalmoscope Basic, <b>MacroView</b> Otoscope Basic, Dispenser, CSM/SPOT 44 Mount - AU Power ( <b>Spot</b> 4400 Monitor not included)
777-PM3WCX-6	Wall System - 77716, <b>PanOptic</b> Plus Ophthalmoscope, <b>MacroView</b> Plus Otoscope, <b>SmartBracket</b> phone adaptor for <b>iExaminer</b> App, Dispenser, CSM/ <b>Spot</b> 44 Mount - AU Power ( <b>Spot</b> 4400 Monitor not included)
BATT22	<b>Connex</b> Spot, 2-Cell LI-Ion Battery
106275	Micro USB Cable
REUSE-11	Cuff, WA, Reusable, Adult
REUSE-12	Cuff, WA, Reusable, Large Adult
4500-35	Fast BP Hose with <b>FlexiPort</b> , 304.8 cm (10 ft)

## DISCOVER THE DIFFERENCE TODAY AT [HILLROM.COM.AU](http://HILLROM.COM.AU)

\*See a full list of parts and accessories online at [hillrom.com.au](http://hillrom.com.au).

### References

1. 2003 std-Alpert BS, Validation of the Welch Allyn Spot Vital Signs blood pressure device according to the ANSI/AAMI SP10: 2002. Accuracy and cost-efficiency successfully combined. *Blood Press Monit.* 2007 Oct;12(5):345-7.
2. Muntner, Paul et al. "Measurement of Blood Pressure in Humans: A Scientific Statement From the American Heart Association." *Hypertension* 73, no. 5. (May 1, 2019). Doi:10.1161/hyp.000000000000087.
3. Meccariello M, Perkins D, Quigley LG, Rock A. Vital Time Savings. *Eval Use an Autom Vital Signs Doc Syst a Medical/Surgical Unit JHIM.* 2010;24(4):46-51.
4. Robert Smith, MD. Blood pressure averaging methodology: Decreasing the risk of misdiagnosing hypertension.

**For more information, please contact your Baxter representative or call Baxter Customer Care on 1800 445 576.**

## **ALWAYS READ THE LABEL AND FOLLOW THE DIRECTIONS FOR USE**

**CONTACT YOUR BAXTER REPRESENTATIVE FOR THE LATEST VERSIONS**

Hillrom reserves the right to make changes without notice in design, specifications and models. The only warranty Hillrom makes is the express written warranty extended on the sale or rental of its products.

Baxter, Connex, Hillrom, iExaminer, Green Series, PanOptic, SureTemp, FlexiPort, SmartBracket, Spot Vital Signs and Welch Allyn are trademarks of Baxter International Inc. or its subsidiaries. Nonin is a trademark of Nonin Medical, Inc.

AUTHORISED AUSTRALIAN SPONSOR Welch Allyn Australia Pty. Limited. 1 Baxter Drive, Old Toongabbie NSW 2146.

[hillrom.com.au](http://hillrom.com.au)

**Baxter**

ANZ-FLC157-230011\_r1 10/2023