



STURDY FUNCTIONALITY WITH ERGONOMIC DESIGN

TECHNICAL SPECIFICATIONS

Height adjustment	730 mm – 1080 mm (without rails)
Sizes	Width: 960 mm Depth: 490 mm
Surface area	Fixed Model Table: 960 mm x 450 mm Tilting Model Table: 740 mm x 450 mm / fixed portion 220 mm x 450 mm
Weight	18.5 kg
Maximum load	16 kg
Material	SURFACE AREA High Pressure Laminate (HPL) UNDERCARRIAGE Metallic grey epoxy coated steel (with plastic cover for H-shaped foot)
Casters	4 x 50 mm double-band with 2 brakes
Model TA270 TT	30° locked positions / 360° full rotational

HPL SURFACE COLOURS

 M91 Hillrom Blue	 M32 Crystal White
 M89 Scandic Beech Dark	 M40 Light Grey
 M47 Steamed Beech	 M95 Light Beige
 M88 Bloomed Cherry	

- The Autolift™ system from Hillrom helps reduce the risk of entrapment and overturning. With a simple push, the surface lifts and moves securely into place.
- The H-shaped base features a reduced width front wheelbase which delivers maximum stability and is compatible in all care environments.
- Locking dual-directional tilted surface provides an ideal reading position.



High Pressure Laminate Surface



To learn more, please contact your Hillrom Representative or Hillrom Customer Care on 1800 445 576.

Always read the labels and instructions for use.

This document is destined solely for use by healthcare professionals. Devices shown in this brochure are intended for use with patients in departments of healthcare establishments. Hill-Rom recommends you carefully read the detailed instructions for safe and proper use included in the documents accompanying the devices. The personnel of the healthcare establishments are responsible for the proper use and maintenance of these devices.

Furniture manufacturer: Hill-Rom SAS, B.P. 14 - Z.I du Talhouët, 56330 Pluvigner - France (also European Representative)

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale of its products.