



## Unique Characteristics

# PROGRESSA® BED SYSTEM

The Progressa bed system responds to the evolving needs of caregivers, patients and health care facilities. It's more than just a bed – this therapeutic device acts as a seamless extension of the health care team. It is designed to support early mobility, and to treat or prevent pulmonary or other complications associated with immobility. It's also flexible enough to configure based on your needs today, and upgrade as those needs change.

### Choose the surface option that's right for your patient:

- **Therapy Surface** (powered air) helps prevent skin ulcers.
- **Pulmonary Surface** (powered air) helps prevent ventilator associated complications with CLRT.

### UNIQUE CHARACTERISTICS:

- StayInPlace® technology helps keep patients from sliding toward the foot of the bed. It also minimises shearing and friction to protect patients' skin.
- Full chair position supports mobility and allows patients to exit out the foot-end of the bed. This full-seated position means you can acclimate immobile patients without moving them out of the bed.
- Cushion-based lateral rotation, percussion and vibration therapy assist in breaking up secretions for immobilised patients. (Cushion features accommodate up to 227 kilos).
- Hands-free foot pedals help you raise and lower the bed when your hands are busy helping your patient.
- Chest X-ray cassette zip supports daily chest X-rays without affecting patients' skin integrity or ripping the surface.
- The second-generation surface uses horizontal air flow that is not impeded by bed articulation (unlike competitors' traditional pinpoint low air loss).
- A built-in oxygen tank holder helps safely transport critically ill patients.
- Bed Exit Alarm.

To learn more, please contact your Hillrom Representative or Hillrom Customer Care on 1800 445 576.

**Always read the labels and instructions for use.** Refer to user manual 1171528 Progressa® Bed Manual.

AUTHORISED AUSTRALIAN SPONSOR Hillrom Pty Ltd. Suite 4.01, 2-4 Lyonpark Road, Macquarie Park 2113 NSW Australia

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale of its products.