



Hillrom™

Welch Allyn
SureTemp® Plus 692
Electronic Thermometer



SureTemp Plus 692 Thermometer

The speed and accuracy you require, and the safety and security you demand.

Easy:

Intuitive user interface, easy-to-read LCD screen and flashing icons, helps select and identify temperature measurements and modes.

Safe:

Designed with interchangeable oral/axillary (blue) and rectal (red) probe wells to address patient safety and help reduce the risk of cross-contamination.

Fast and Accurate:

With proven accuracy, the SureTemp 692 thermometer captures quick readings—approximately 4-6 second oral, 10-12 second paediatric axillary and rectal and 13-15 second adult axillary measurements.

Secure:

To top it all off, our security features make sure the unit stays where it is supposed to be. A location ID on the display and other innovative features help prevent theft.

- Easy-to-use, intuitive operation
- Interchangeable oral/axillary and rectal probe wells help reduce the risk of cross-contamination

SureTemp Plus 692 Thermometer

Setting the standard of care



The SureTemp Plus 692 with standard wall mount has storage for an additional 50 probe covers.

Interchangeable rectal probe and well reduce risk of cross-contamination

The SureTemp Plus 692 Thermometer

Setting the standard of care even higher

SureTemp Plus 692 Product Specifications

Calibration Accuracy	±0.1 °C Monitor Mode	
Temperature Predict Time	Oral: Adult Axillary (18 years and older): Paediatric Axillary (17 years and younger): Rectal:	Approximately 4-6 secs Approximately 15 secs Approximately 10 secs Approximately 10 secs
Unit Dimensions	215 mm x 81 mm x 62 mm	
Weight	357 g	
Ambient Temperature Range	10 °C to 40.0 °C	
Patient Temperature Range	26.7 °C to 43.3 °C	
Power Source	Ⓢ 1.5Vdc AA batteries	
Battery Life	Approximately 5,000 temperatures	
Warranty Period	3 Years	

Ordering Information:

THERMOMETERS

01692-400	SureTemp Plus 692 Electronic Thermometer, Wall Mount, Security System with ID Location Field, Oral Probe and Oral Probe Well with 1.2 m cord
01692-201	SureTemp Plus 692 Electronic Thermometer, Wall Mount, Security System with ID Location Field, Rectal Probe and Rectal Probe Well with 1.2 m cord
01692-500	SureTemp Plus 692 Electronic Thermometer, Wall Mount, Security System with ID Location Field, Oral Probe and Oral Probe Well with 2.7 m cord

ACCESSORIES

21326-0000	Wall Holder
21330-0000	Holder for Spare Probe and Well
06138-000	Calibration Key Assembly, 690/692
02892-000	Probe and Well Kit, 1.2 m Rectal
02892-100	Probe and Well Kit, 2.7 m Rectal
02893-000	Probe and Well Kit, 1.2 m Oral
02893-100	Probe and Well Kit, 2.7 m Oral
02891-0000	Probe Well Blue
02891-1000	Probe Well Red
08273-000	Thermometer Stand with Holder for Spare Probe and Well
21333-0000	Bracket, M600 Stand Upgrade, with Holder for Spare Probe and Well

REPLACEMENT PROBE COVERS

05031-101	Case of 1,000 Probe Covers
05031-105	Case of 5,000 Probe Covers



Hillrom™

To find out more, please contact your Hillrom Sales Representative or Customer Service on 1800 445 576.

Always read the labels and instructions for use. Contact your Hillrom Representative for the latest version.