



Hillrom™

# Welch Allyn® Multi-Parameter Vital Signs Devices



Spot Vital Signs® 4400 Device



Connex® Spot® Monitor



Connex Vital Signs Monitor



Connex Integrated Wall System

Workflows	Spot	Spot, Office, Intervals	Spot, Office, Intervals, Continuous	Spot, Office, Intervals, Continuous
Device Type	Basic	Intermediate	Advanced	Advanced
Patient Populations	Adult, Pediatric	Adult, Pediatric, Neonatal*	Adult, Pediatric, Neonatal*	Adult, Pediatric, Neonatal*
Device Alarms	—	Yes	Yes	Yes
Third-Party Alarm Notification	—	—	Yes	Yes
Parameters	Blood Pressure, Temperature, SpO <sub>2</sub>	Blood Pressure, Temperature, SpO <sub>2</sub> , RRp	Blood Pressure, Temperature, SpO <sub>2</sub> , etCO <sub>2</sub> , SpHb, ECG, RRa	Blood Pressure, Temperature, SpO <sub>2</sub> , SpHb, ECG, Otoscopy, Ophthalmoscopy
SureBP® Technology (15-sec)	Yes	Yes	Yes	Yes
Blood Pressure Averaging	Yes	Yes	Yes	Yes
Diagnostic Handles	—	—	—	Otoscope, Ophthalmoscope
SpO <sub>2</sub> Options	Nonin®	Nellcor®, Masimo®, Nonin	Nellcor, Masimo	Nellcor, Masimo
Temperature Options	SureTemp® Technology	SureTemp Plus, Braun PRO 6000	SureTemp Plus, Braun PRO 6000	SureTemp Plus, Braun PRO 6000
Early Warning Scoring	—	Up to 6 protocols	Up to 6 protocols	Up to 6 protocols
Custom Documentation Fields	5 Fixed and Selectable	Up to 20	Up to 20	Up to 20
Wireless Connectivity	—	WiFi®, Bluetooth®	WiFi	WiFi
Wired Connectivity	USB	USB, Ethernet	USB, Ethernet	USB, Ethernet
Display	7" Color LCD, Touchscreen	7" Color LCD, Touchscreen	8.9" Color LCD, Touchscreen	8.9" Color LCD, Touchscreen
Scale Integration	—	Yes	Yes	Yes
Mounting Options	Desktop Stand, Classic Stand	Power Management Stand, Classic Stand, Work Surface Stand, Wall Mount, Integrated Wall System, Desktop Stand	Cable Management Stand, Classic Stand, Wall Mount	Wall Mounted Unit
Accessory Options	—	Scales, Barcode, RFID	Scales, Barcode, RFID	Scales, Barcode, RFID
Memory (# Readings)	50	400	400	400

\*Not all parameters may be used with neonates, plus refer to appropriate IFU for details.

For more information, please contact your Hillrom representative at 1-800-535-6663.

[hillrom.com](http://hillrom.com)

Hillrom reserves the right to make changes without notice in design, specifications and models. The only warranty Hillrom makes is the express written warranty extended on the sale or rental of its products.