



CONNEX® VITAL SIGNS MONITOR

6	X	X	X	X	X
MODEL	RADIO	PULSE OXIMETRY	RESPIRATION	THERMOMETRY	PRINTER
6	7 = Not Included	N = Nellcor® SpO ₂	C = CO ₂	T = SureTemp® Plus	P = Internal
	8 = Included	M = Masimo® SpO ₂	R = RRa*	E = Braun® PRO 6000	X = Not included
			X = Not included	D = Both	
				X = None	

* RRa available as upgrade only.



CONNEX® INTEGRATED WALL SYSTEM

8	X	X	X	X	X	X
MODEL	RADIO	PULSE OXIMETRY	ORAL TEMP	STORAGE	EAR TEMP	HEADS
8	4 = Not Included	N = Nellcor® SpO ₂	T = SureTemp Plus	V = SCS (standard cuff storage)	E = Braun® PRO 6000	C = MacroView® + Coaxial Ophthalmoscope
	5 = Included	M = Masimo® SpO ₂	X = None		X = None	P = MacroView + PanOptic 118
		X = None				2 = MacroView + PanOptic Basic
						3 = MacroView + PanOptic Plus
						X = None

POWER CORDS

For Connex devices, power cords are indicated as “-X” at the end of the part number as a 1-digit alphanumeric code (e.g., 84MTVEC-B or 67MXTP-B).

POWER CORD KEY	
B	North America
2	Europe
4	United Kingdom
6	Australia/New Zealand
7	South Africa

MODEL	67	68	84	85
WIRELESS CONNECTIVITY				
Upgradeable	✓		✓	
Integrated Radio		✓		✓
WIRED CONNECTIVITY				
USB	✓	✓	✓	✓
Ethernet	✓	✓	✓	✓
Nurse Call	✓	✓	✓	✓
NUMBER OF ACCESSORIES CAPABLE OF RUNNING	4	4	4	4

For more information, please contact your local distributor or Hillrom sales representative at 1-800-535-6663.

hillrom.com

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.